



# Frequently asked questions about Change Healthcare Cardiology Hemodynamics™

System capabilities, outputs, and configurations

Our integrated monitoring system, Change Healthcare Cardiology Hemodynamics, aggregates haemodynamic data, waveforms, and images into one cardiac patient record, eliminating redundant data entry and improving care decisions.

Continuous vitals monitoring provides critical information during Cath procedures.

**Optum**

**Q: What are the capabilities of the solution?**

A: Change Healthcare Cardiology Hemodynamics is an integrated haemodynamic monitoring system for monitoring vital signs and performing measurements and calculations. The system documents procedure and patient data, and interfaces with other systems and devices before, during, and after procedures.

The system provides comprehensive haemodynamic monitoring and recording for cardiology, cardiac catheterisation, electrophysiology, radiology, and invasive radiology procedures.

Change Healthcare Cardiology Hemodynamics also imports patient data from other hospital information systems and shares haemodynamic data in return. Change Healthcare Cardiology Hemodynamics seamlessly interfaces with hospital information systems as well as image management, archiving, and reporting systems via Change Healthcare Cardiology solution.

**Q: What patient vitals data does the system monitor?**

A: Change Healthcare Cardiology Hemodynamics is intended for physiological/haemodynamic monitoring, as well as clinical data acquisition. Medical image and data processing, and analytical assessment is done via the Change Healthcare Cardiology solution.

Change Healthcare Cardiology Hemodynamics provides patient monitoring via:

- ECG leads for ECG traces, heart rate, and respiration rate
- Invasive Blood Pressure (IBP) using up to 4 simultaneous IBP channels, plus optional Fractional Flow Reserve (FFR) measurement
- SpO2 finger clip for peripheral saturation and Heart rate (HR)
- Non-invasive blood pressure (NIBP) cuff
- Body temperature probe
- Thermal Dilution Cardiac output temperature probe
- Capnograph and etCO2 for inspired CO2 and respiration rate

**Q: What are the components of the system?**

A: Change Healthcare Cardiology Hemodynamics is composed of:


- **A control and documentation unit** that performs measurements, records full disclosure, and takes samples. Procedure notes are entered directly into this unit, and overall data input and management of the patient and procedure data is performed within this unit.
- **The ‘front-end’ unit and the real-time monitor.** The clinical unit acquires, analyses, and displays patient vital signs and other pertinent clinical data. The clinical data is subsequently displayed on monitors..

**Q: How is the system set up and integrated?**

A: Providers can use the Change Healthcare Cardiology Hemodynamics system in one of three configurations.

## Configuration 1

Hemo + Change Healthcare Cardiology CVIS (Cath and EP reports) + CPACS




This configuration provides a comprehensive solution enabling:

- Extensive interface (import and export) capabilities with the EPR/PAS, modalities, and various consumers (e.g., ADT, orders, lab results, allergies, medications, and reports)
- The sharing of patient and procedural data between multiple Haemodynamics stations
- Monitoring and documenting invasive procedures via Change Healthcare Cardiology Hemodynamics
- Storing, managing, and viewing images captured during the procedure via the different modalities, for example, by Cath and IVUS
- Producing physician reports

## Configuration 2

Network Hemo + Change Healthcare Cardiology CPACS



This configuration provides a haemodynamic monitoring solution, which enables:

- Extensive interface (import and export) capabilities with the EMR, modalities, and various consumers (e.g., ADT, orders, lab results, allergies, medications, procedure results, and reports)
- The sharing of patient and procedural data between multiple Haemodynamics stations
- Monitoring and documenting invasive procedures via Change Healthcare Cardiology Hemodynamics

## Configuration 3

Stand-alone Hemo



This configuration focuses solely on the haemodynamic monitoring solution, which enables a very limited interface and file export.

## Q: What is the cardiology procedure workflow using this system?

1. The technologist places a generic EP or Cath Lab CV order in the EMR scheduling system.
2. Change Healthcare Cardiology Hemodynamics will refresh the worklist to display this order.
3. Cath, EP (and/or other modalities) will query the DICOM Worklist for ordered procedure.
4. Patient and procedure data is automatically enhanced with additional interfaced information, such as lab results and EPR/PAS data **(Configurations 1 and 2)**.
5. The technologist will launch Change Healthcare Cardiology Hemodynamics for the relevant procedure and begin monitoring and documenting the procedure.
6. Procedural data is populated to the physician report **(Configuration 1)**.
7. The technologist will end the procedure on the Cath Lab modality (or any other modalities, such as ultrasound or intravascular ultrasound).
8. Procedure data from the modality will be stored to Change Healthcare Cardiology CPACS **(Configuration 1)**.
9. The EP modality will send the XML file to Change Healthcare Cardiology CVIS **(Configuration 1)**.
10. Charges and the documentation of used devices will be verified by the technologist in Change Healthcare Cardiology Hemodynamics **(Configurations 1 and 2)**.
11. The technologist will sign the Hemo procedure in Change Healthcare Cardiology Hemodynamics and generate a PDF.
12. As all cardiovascular images are now accessible, the cardiologist/electrophysiologist will launch the physician report and viewer to view images and complete the report. Once signed, a structured report PDF is produced in Change Healthcare Cardiology **(Configuration 1)**.
13. Additional interfaces are triggered to export procedural data to the EPR/PAS system in HL7 or XML format **(Configurations 1 and 2)**.



[Visit us to learn more](#)

# Optum

optum.com

Optum is a registered trademark of Optum, Inc. in the U.S. and other jurisdictions. All other brand or product names are the property of their respective owners. Because we are continuously improving our products and services, Optum reserves the right to change specifications without prior notice. Optum is an equal opportunity employer.

© 2023 Optum, Inc. All rights reserved. WF11302761 09/23



Change Healthcare Canada Company  
10711 Cambie Road  
Richmond, BC V6X 3G5  
Canada

CPR-0020084 Rev 2