



# Change Healthcare Cardiology

## ECG Management™

### Vendor Compatibility List

Product Documentation

**Vendor Compatibility CHC ECGM**

Document ID: CGUD-10113

Revision 1.0

Effective Date: January 2021

**Rx Only**

# Document and release information

## Copyright and trademarks

Copyright 2021 Change Healthcare LLC and/or one of its subsidiaries. Change Healthcare and the Change Healthcare logos, Change Healthcare Cardiology, Change Healthcare Cardiology Hemo, Change Healthcare Cardiology ECG Management, Change Healthcare Cardiology Customizer, Change Healthcare Management Console, Change Healthcare Cardiology Cath, Change Healthcare Cardiology Reference Viewer, Change Healthcare Cardiology Charge Manager, Change Healthcare Cardiology Echo, Change Healthcare Cardiology EP, Change Healthcare Cardiology Web, Change Healthcare Cardiology Holding Area Charting, Change Healthcare Cardiology Inventory Admin, Change Healthcare Cardiology Nuclear Medicine, Change Healthcare Cardiology Nurse Charting Solution, Change Healthcare Cardiology Reporting, Change Healthcare Cardiology Analytics and Change Healthcare Cardiology Statistical Report Center are trademarks of Change Healthcare LLC. All other product or company names mentioned may be trademarks, service marks or registered trademarks of their respective companies. All Rights Reserved.

While reference may be made throughout to other entities, all software and documentation for this and previous versions of this software are owned by Change Healthcare and/or one of its subsidiaries.

## Usage

Use of this documentation and related software is governed by a license agreement. This documentation and related software contain confidential, proprietary, and trade secret information of Change Healthcare LLC and/or one of its subsidiaries or affiliates and is protected under United States and international copyright and other intellectual property laws. Use, disclosure, reproduction, modification, distribution, or storage in a retrieval system in any form or by any means is prohibited without the prior express written permission of Change Healthcare LLC. This documentation and related software are subject to change without notice.

All patient information used in this manual is fictitious. The information does not refer to any persons either living or deceased.

Change Healthcare LLC and/or one of its subsidiaries (Change Healthcare) assumes no responsibility for the use or suitability of software or other equipment not provided by Change Healthcare or its affiliated companies. Although Change Healthcare strives to ensure the accuracy of the provided instructions, neither Change Healthcare nor its affiliated companies accept any liability or damages for omissions contained in this manual. This document does not imply a legal contract. Readers assume responsibility for the use and interpretation of the information contained herein.

## Product

Change Healthcare Cardiology ECG Management

## Publication number and revision history

Document ID: CGUD-10113


Revision	Date	Description
1.0	January 2021	Initial version

## Corporate address

Change Healthcare Holdings LLC  
100-3055 Lebanon Pike  
Nashville, TN 37214  
U.S.



Change Healthcare Canada Company  
10711 Cambie Road  
Richmond, BC, V6X 3G5  
Canada



Authorized Representative in the European Community  
Emergo Europe B.V.  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

## Initial Importer/Distributor for the United States

Change Healthcare Technologies LLC  
220 Davidson Avenue  
Somerset, NJ 08873-4149  
U.S.

## Importer for Europe

Change Healthcare Ireland Solution Limited  
1 Woodford Business Park  
Tumapin Lane  
Santry, Dublin 17, Ireland

## Authorized regulatory sponsor in Australia

Emergo Australia  
201 Sussex Street  
Darling Park, Tower II  
Level 20 Sydney, NSW, 2000  
Australia

**Authorized regulatory sponsor in New Zealand**

CARSL Consulting  
50B Puketapu Road  
Taradale  
Napier 4112  
New Zealand

**Conformity****Support information**

For technical support or any service related to your product, call the toll-free number for your region listed at <https://iwcrm.changehealthcare.com>. Standard coverage may include software support, hardware support, and software updates as covered in the support maintenance agreement.

**Reader comments**

Comments or suggestions regarding this publication are welcome. Contact Change Healthcare Support at <https://iwcrm.changehealthcare.com/>.

# Table of contents

- About this document..... 1
- Intended audience..... 1
- 1 Chapter 1 - Rest ECG..... 2
  - 1.1 Datamed™ ..... 2
    - 1.1.1 DatamedFT™ ..... 2
    - 1.1.2 DataMedWL™ ..... 5
  - 1.2 Mortara ELI™ ..... 5
- 2 Chapter 2 - Stress ECG ..... 6
- 3 Chapter 3 - Holter ECG carts ..... 7

# About this document

This document provides the list of ECG Modality vendors to which Change Healthcare Cardiology ECG Management is compatible.

This document is applicable to Change Healthcare Cardiology 14.1.1, 14.2 and 14.3.

## Intended audience

This document is written for Change Healthcare Cardiology ECG Management users who wish to know to which ECG modalities are compatible.

Typically, the users are:

- Physicians
- Administrators
- Technicians
- Other medical personnel

# 1 Chapter 1 - Rest ECG

## 1.1 Datamed™

### 1.1.1 DatamedFT™

Supported modalities:

Vendor	Supported Devices
AMEDTEC Medizintechnik Aue GmbH	AMEDTEC ECGPro
Cardiolex	EC Sense
Cardioline SpA	Cardioline® Rest ECG cardiographs (including ClickECG)
corpuls® GS	Elektromedizinische Geräte G. Stemple GmbH corpuls3 defibrillator/monitors
custo med® GmbH	custo diagnostic, cardio 100, cardio110, cardio 100BT, cardio 110BT, cardio 130, cardio 200BT
Diagnostica AS	MeditConcept cardiograph
Dräger	Infinity® patient monitors (Delta, Delta XL, Kappa, Omega Package, Omega-S Package, Acute Care System, M540)
EB Neuro SpA (Ates Medica Device)	Easy ECG series cardiographs
Edan Instruments, Inc.	SE-12, SE-12 Express, SE-1200, SE-1200 Express, SE-1201, SE-15, SE-1515, SE-16, SE-18, SE-3 series, SE-300, SE-301, SE-601 series, and SE-1010 cardiographs
Esaote®	P80, P80 Power, P8000, P8000 Power cardiographs AT-10 plus cardiographs MyTrace 16 cardiographs MyCardioPad and MyCardioPad XL cardiographs
Fukuda Denshi	FCP-7101, FX-7102, FX-7202, FX-7302, FX-7542, FX-8222, and FX-8322 cardiographs FX-7402 cardiographs Patient monitors with the CVW Monitoring Gateway
GE® Healthcare	GE® MAC®-600 cardiographs GE® MAC®-800, MAC®-1600, MAC®-3500, MAC®-5500, and MAC®-5500HD cardiographs GE® MAC®-1200 cardiographs GE® MAC®-2000 cardiographs

Vendor	Supported Devices
	GE® MAC®-5000 cardiographs GE® CardioSoft™ cardiographs GE® DASH®, SOLAR®, and CARESCAPE® patient monitors Marquette® MAC®-VU cardiograph
<i>Gem-Med, s.l.</i>	Gem Heart cardiographs
<i>Hillrom (including Cardiac Science, Mortara®, and Welch Allyn®)</i>	Mortara® ELI® series cardiographs, Welch Allyn® Connex Cardio™ Welch Allyn® CP150 with DICOM® Option Welch Allyn® CP-series cardiographs with CardioPerfect® Workstation Quinton® Eclipse™ 850, Eclipse™ Plus, and Eclipse™ Premier cardiographs Burdick® Atria 3000, Atria 3100, Atria 6100, Atria 8300, and Atria 8500 cardiographs Burdick® CareCenter MD™ cardiographs
<i>HP (Philips)</i>	PageWriter® M1700A (XLi), M1701A (XLs), M1702A (XLe) cardiographs PageWriter® M1770A (100/100i), M1771A (200/200i), M1772A (300/300i/300pi) cardiographs PageWriter® TC70 Card Guard CG-7000DX BT
<i>Mindray®</i>	BeneHeart® R3 and R12 cardiographs BeneHeart® D3 and D6 defibrillator/monitors V-Series monitors Passport 2 and Spectrum monitors with View 12 ECG Module connected to Panorama Central Patient monitors and defibrillators with 12 ECG Module connected to Hypervisor Central Station (v6.6 or higher) Passport 12m and 17m monitors with 12 ECG Module connected to DPM Central Station (v01.01 or higher) Passport 12m and 17m monitors with 12 ECG Module connected to BeneVision® Central Station (v02.01.01 or higher)
<i>Nihon Kohden® Corporation</i>	Cardiofax® ECG-1250A/K, ECG-1350A/K, ECG-1550A/K, ECG-2350, and ECG-2550 Cardiofax® ECG-1150A/K, ECG-9000 series (9020, 9130K, 9130P, 9230K, 9232K, 9320, etc.) Polaris One system Lifescope® BSM-series patient monitors with QP-992P Gateway Lifescope® BSM-series patient monitors with QP-984PK Gateway



<b>Vendor</b>	<b>Supported Devices</b>
<i>Norav Medical™ Ltd</i>	PC ECG 1200™ family (1200M™, 1200HR™, 1200S™, 1200W™, Blue-ECG™), Green ECG LAN™, NHH 1200™, NECG-3™, NECG-12™, NECG-12C™
<i>Ortivus AB</i>	Mobimed cardiographs, Mobimed Smart, CoroNet patient monitors
<i>Philips® Healthcare</i>	PageWriter® Touch, PageWriter® Trim II, PageWriter® Trim III, and PageWriter® TC (10/20/30/50/70) cardiographs  IntelliVue® (MP5/MP30/MP50/MP70/MP90/X2) patient monitors with 12-Lead option  HeartStart® MRx defibrillator/monitors  LIFEPAK® 12 and LIFEPAK® 15 defibrillator/monitors
<i>QRS Diagnostic</i>	Office Medic™
<i>SCHILLER AG</i>	CARDIOVIT® AT-102 G2 and CS-104 cardiographs  CARDIOVIT® AT-101 series and AT-102 cardiographs  CARDIOVIT® AT-102 Plus, MS-2007, MS-2010, MS-2015, AT-170 cardiographs  CARDIOVIT® CS-200, CS-200 Office, CS-200 Touch, and CS-200 Excellence cardiographs  DEFIGARD® DG-5000, ARGUS PRO LifeCare, and ARGUS PRO LifeCare2 with ARGUS PRO  CARDIOVIT® FT-1  SEMA ECG Management System  SEMA3 ECG Management System  CARDIOVIT® CS-6/12, AT-3, AT-4, AT-5, AT-6, AT-60, CS-100, AT-10, and AT-2 plus cardiographs  CARDIOVIT® AT-104 PC, and CardioLaptop AT-110 cardiographs  CARDIOVIT® AT-10 plus cardiographs
<i>SEIVA s.r.o.</i>	ECG Praktik, ECG Praktik Portable, ECG Cardio Writer, ECG Cardio Touch
<i>Spacelabs Healthcare</i>	CardioExpress® cardiographs  XPREZZON®, Qube®, Ultraview, and UltraviewSL® patient monitors with ICS and ICS G2
<i>VectraCor®, Inc.</i>	VectraplexECG™ System
<i>Verathon® Inc</i>	Heartscape™ 3D ECG System, Prime™ ECG cardiographs
<i>ZOLL® Medical Corporation</i>	E-Series® (prior to rev 4.01) and M-Series® defibrillator/monitors  X-Series® defibrillator/monitors

### 1.1.2 DataMedWL™

Supported modalities:

Vendor	Supported Devices
<i>Edan Instruments, Inc.</i>	SE-1201, SE-1200E(Basic), SE-1200E, SE-12E, SE-15, SE-1515, SE-16, and SE-18 cardiographs
<i>GE® Healthcare</i>	MAC® 800, MAC® 2000, MAC® 3500, MAC® 5000, MAC® 5500, and MAC® 5500HD
<i>Hillrom (including Cardiac Science, Mortara®, and Welch Allyn®)</i>	ELI® series cardiographs Welch Allyn® CP150 with DICOM option Burdick® Atria 3100/6100 and 8300/8500 cardiographs Eclipse™ Premier cardiographs
<i>Nihon Kohden® Corporation Devices: Cardiofax®</i>	ECG-1250A/K, ECG-1350A/K, ECG-1550A/K, ECG-2350, and ECG-2550 cardiographs
<i>Philips® Healthcare</i>	PageWriter® Touch, PageWriter® TC cardiographs
<i>SCHILLER AG</i>	CARDIOVIT® AT-102, AT-102 Plus, AT-10 plus, AT-170, FT-1, AT-102 G2, MS-2007, MS2010, and MS-2015 cardiographs CARDIOVIT® CS-104, CS-200, CS-200 Office, CS-200 Touch, and CS-200 Excellence cardiographs

### 1.2 Mortara ELI™

The supported modalities include the following devices: ELI10®, ELI150®, ELI250®, ELI350® and ELI380®.

---

**Note:** The Mortara® product line is now part of the Welch Allyn® cardiology portfolio.

---

## 2 Chapter 2 - Stress ECG

Supported modalities:

- Mortara X-Scribe®
- GE Case®
- HP – Philips ST80i®
- Quinton – Q-Stress®

---

**Note:** ECG Management supports the transfer of DICOM encapsulated PDFs. DICOM import is supported only if the DICOM file adheres to the standard DICOM encapsulated PDF format, and the DICOM contains a tag indicating the file as a Stress procedure.

---

---

**Note:** The Mortara® product line is now part of the Welch Allyn® cardiology portfolio.

---

## 3 Chapter 3 - Holter ECG carts

The supported modality is: Mortara H-Scribe®

---

**Note:** ECG Management supports the transfer of DICOM encapsulated PDFs. DICOM import is supported only if the DICOM file adheres to the standard DICOM encapsulated PDF format and the DICOM contains a tag indicating the file as a Holter procedure.

---

---

**Note:** The Mortara® product line is now part of the Welch Allyn® cardiology portfolio.

---